



Quick reference guide

Infinity[®] Acute Care System

Infinity M540 patient monitor
Software VG7.n

This Quick reference guide is not a replacement or substitute for the Instructions for Use.

Any use of the medical device requires full understanding and strict observation of the Instructions for use.

Contents

Contents	2
Trademarks	3
Side view of the M540	3
Dialogs	4
Calibrating the touchscreen	4
To select a view	4
Discharging a patient	4
Selecting a new patient category	5
To adjust the alarm volume	5
Pausing acoustic alarm signals (audio pause)	5
Pausing alarms at the M540	6
Pausing alarm monitoring temporarily	7
Activating or deactivating alarm monitoring	7
Setting the upper and lower alarm limits	8
Archive function	8
To access stored events	9
Electrode placement	10
ECG parameter setup functions	12
Respiration parameter setup functions	14
SpO ₂ and Pulse CO-Ox monitoring with Masimo SET MCable	16
SpO ₂ and pulse rate monitoring with Nellcor OxiMax MCable	17
Temperature parameter setup functions	18
Non-invasive blood pressure parameter setup functions	19
To assign a pressure label manually	20
Zeroing a specific transducer	20
Invasive pressure parameter setup functions	21
Mainstream CO ₂ monitoring	22
Microstream CO ₂ monitoring	23

Trademarks

- Infinity®
- MCable®
- MPod®

are trademarks of Dräger.

- Masimo®
- SET®
- Pulse CO-Oximeter Signal Extraction®

are trademarks of Masimo corporation.


- Microstream®

is a trademark of Oridion Medical 1987 Ltd.

- Nellcor®
- OxiMax®
- SatSeconds®

are trademarks of a Medtronic company.

Side view of the M540



- A Non-invasive blood pressure connector
- B Temp (2) / Aux connector
- C **SpO₂** connector
- D Hemo connector
- E CO₂ connector
- F **ECG** connector



Dialogs

To access parameter-specific setup pages directly, touch the corresponding parameter fields on the main screen.

Calibrating the touchscreen

If the touchscreen is out of alignment, the user can calibrate it at any time.

To calibrate the touchscreen

- 1 Touch the **Menu** function key.
- 2 Touch the **Screen setup** tab > **Settings** tab > **Touch calib.**
- 3 Touch each cross appearing successively in each corner of the screen.
or
- 1 Push and hold the following two keys simultaneously:  
- 2 Touch each cross appearing successively in each corner of the screen.

To select a view

- Touch the currently selected function key several times (for example, **View 5**) to scroll through the available view labels.

Discharging a patient

A patient discharge has the following effect at the M540:

- All demographic data are removed from the screen
- Any active recordings are canceled at the Cockpit if the M540 is docked in an IACS configuration
- Factory or user default limit settings are restored
- The message **Discharged, Touch screen to resume monitoring** appears

To discharge a patient

Press the **Discharge** function key (if available for display).

or

- 1 Press the **Menu** function key.
- 2 Touch the **Main** tab.
- 3 Touch **Discharge**. The message **Caution discharge will delete patient data** appears.
- 4 Touch **Discharge**.

Selecting a new patient category

After selecting the patient category, the new patient category label appears in the header bar.

A patient category change does not affect the following settings: the patient and physician names, patient ID, birth date, admit date, height, and weight.

To select a new patient category

- 1 Press the **Menu** function key.
- 2 Touch the **Patient setup** tab.
- 3 Touch **Patient category** and then select the appropriate category (**Adult**, **Pediatric**, **Neonate**). The message **Changing category will change alarm settings and algorithmic processing** appears.
- 4 Touch **OK**.

To adjust the alarm volume

- 1 Press the **Alarms** function key.
- 2 Touch the **Alarm volume** button and select the desired volume level (**Off**, 5%, 10 to 100% in increments of 10%).

When the M540 is on wireless transport in an IACS configuration or in standalone mode and it loses its connection to the ICS, the alarm volume setting is automatically set to 100%. The setting **Off** is no longer available until the connection to the ICS is restored.

Pausing acoustic alarm signals (audio pause)


Active alarms can be paused, or silenced, at the M540 for two minutes. In addition to pausing alarms, the setting of the quiet mode feature determines how subsequent alarm conditions are announced.

You can initiate an audio pause in several ways:

- From an M540 in standalone mode or on wireless transport
- From an ICS
- From the remote view of another Infinity monitor within the same monitoring unit
- From the Cockpit when the M540 is docked in an IACS configuration
- From a remote device when remote control and remote silence are activated at the remote device and the Cockpit (refer to the IACS Monitoring Applications IFU)
- When first turning on the device

Pausing alarms at the M540

The following happens at the M540 when you pause active alarms:

- All acoustic alarm signals are paused for a maximum of about two minutes.
- The **Audio paused** message appears in the alarm message header of the Cockpit along with the timer and the following symbol: 
- The alarm message appears in the color corresponding to the alarm priority.
- The parameter field no longer flashes in the color corresponding to the alarm priority. It appears in solid color.
- The alarm bar no longer flashes for high-priority and medium-priority alarm conditions.

The behavior of new alarm condition while the systems is in an audio pause state is determined by the **Quiet mode** setting.

To initiate an audio pause at the M540

- Press the yellow key on the M540:



Pressing the key again cancels the audio pause period and all alarm events are reported as usual.

To initiate an audio pause remotely

- Press the following key on the main menu bar of the ICS to audio pause alarms for all assigned patients:



Press the same button in the viewport area to pause alarm tones for an individual patient. For more information, refer to the instructions for use of the ICS.

- Refer to the instructions for use of any remote device within the same monitoring unit for instructions on how to initiate an audio pause.

Pressing the key again cancels the audio pause period and all alarm events are reported as usual.

Pausing alarm monitoring temporarily

If the password-protected alarm pause feature is activated, you can pause alarm monitoring temporarily. The alarm pause duration is adjustable from 1 minute to 5 minutes.

The following happens when you pause alarm monitoring:

- Acoustic and optical alarm signals for new alarm conditions are suppressed for all parameters until alarm monitoring begins again
- Alarm signals for any active alarm condition stop immediately
- The alarming parameter field and alarm bar return to the pre-alarm state
- Alarm messages are removed from the alarm message field in the header bar
- The far right field of the header bar turns yellow and displays the alarm message **All alarms pause**, a timer, and the following symbol:



To pause alarm monitoring temporarily

- 1 Press the **Alarms** function key.
- 2 Touch **All alarms pause**.

As soon as the alarm pause period ends, the M540 generates acoustic and optical alarm signals as needed.

To activate alarm monitoring after pausing

- 1 Press the **Alarms** function key.
- 2 Touch **All alarms pause** before the alarm pause period ends to cancel the alarm pause.

Activating or deactivating alarm monitoring

If the password-protected alarm pause feature is set to **No timeout**, the following happens when you deactivated alarm monitoring:

- All acoustic and optical alarm signals for new alarm conditions are suppressed for all parameters until alarm monitoring is manually activated again
- Acoustic alarm signals for any active alarm condition stop immediately
- The alarming parameter field and alarm bar return to the pre-alarm state
- Alarm messages are removed from the alarm message field of the header bar
- The far right field of the header bar turns yellow and displays a message that all alarms are off and the following symbol appears:



To deactivate alarm monitoring

- 1 Press the **Alarms** function key.
- 2 Touch **All alarms off**.

To activate alarm monitoring after deactivating

- 1 Press the **Alarms** function key.
- 2 Touch **All alarms off**.

The M540 provides acoustic and optical alarm signals again when it detects a new alarm condition.

Setting the upper and lower alarm limits

You can configure the upper and lower alarm limits of a parameter manually to trigger acoustic and optical alarm signals if a parameter goes above or below the set limits.

To set an individual parameter's alarm settings

- 1 Touch the parameter field (for example, **HR**) to access that parameter's dialog.
- 2 Touch the tab for configuring the parameter limits (for example, **HR limits**).
- 3 Touch the upper or lower alarm limit.
- 4 Touch the up or down arrow to change the alarm limit setting.
- 5 Touch **OK**.

Archive function

The archive function setting determines what happens in response to an alarm limit violation. The available settings are:

- **Off** – no event is stored and no recording is generated.
- **Store** – stores the event for later review.
- **Record** – stores the event for later review when in standalone mode or on wireless transport. Once the M540 is docked, any stored event is transferred to the alarm history of the IACS where you can request a manual recording of the event.
- **Str/Rec** – generates a timed recording for an M540 docked in an IACS configuration and stores the event.


When the M540 is in standalone mode or on wireless transport, this setting stores an event for later review. Once the M540 is docked, any stored event is transferred to the alarm history of the IACS where you can request a manual recording of the event.

To configure an individual parameter's archive function

- 1 Touch the parameter field to access that parameter's dialog (for example, **HR**).
- 2 Touch the tab for configuring the parameter limits (for example, **HR limits**).
- 3 Touch **Archive** and toggle to one of the following settings: **Off**, **Store**, **Record**, **Str/Rec**.
- 4 Touch **X** to close the dialog.

To access stored events


Select the **Review** function key. The following diagram shows the **Event recall** dialog.




- A Event list with date, time, and cause of each event.
- B Arrow buttons for scrolling through the event list.
- C **View** button for viewing a single event in greater detail.
- D **Delete** button for deleting an event.
- E **Lock** button for locking an event (if the M540 is docked in an IACS configuration, pressing this button has no effect on events at the Cockpit).

Electrode placement for adult and pediatric patients

Standard configuration, three electrodes (IEC/AHA)




314




315

Standard configuration, five electrodes (IEC/AHA)




316




317

Pacer configuration, five electrodes (IEC/AHA)





318




319


Standard configuration, six electrodes (IEC/AHA)





12-lead configuration, ten electrodes for 12-lead Rest ECG monitoring (AHA)



12-lead configuration, ten electrodes for 12-lead Rest ECG monitoring (IEC)



Electrode placement for neonates



ECG parameter setup functions

All ECG parameter setup functions take place in the **ECG** dialog.

Selection	Available settings	Description
Settings – ECG 1 page		
Tone volume ¹⁾	Off , 5, 10 (default) to 100% in increments of 10%	Selects the volume of the pulse tone. If you dock the M540 in an IACS configuration, this setting is replaced by the pulse tone volume setting of the Cockpit. When you undock the M540, this setting is replaced by the Transport pulse tone setting configured under the Alarm setup tab.
Tone source ¹⁾	ECG (default), PI	Selects the source of the pulse tone.
ECG filter ¹⁾	<ul style="list-style-type: none"> – Off – provides the greatest sensitivity to noise or artifact (the message Filter off appears in the waveform channel) Passband: 0.08 - 40Hz – Monitor (default) – recommended for standard monitoring; reduces wandering isoelectric line, muscle artifact, and power line interference. No message appears in the waveform channel. Passband: 0.5 - 40Hz – ESU – reduces signal distortion during electrosurgery (the message Filter ESU appears in the waveform channel). This selection is not available when the M540 is in standalone mode. Passband: 0.5 - 16Hz <p>12-lead monitoring is not available when the ESU filter is enabled. Likewise, the ESU filter selection is not available when you are using 12-lead monitoring.</p>	<p>Controls the sensitivity to various artifact sources.</p> <p>When the M540 is set to OR alarms and the filter selection is set to Monitoring:</p> <ul style="list-style-type: none"> – the hardware low pass ESU filter is activated. – RRi is unavailable – 12-lead ECG monitoring is unavailable <p>None of these settings are of diagnostic quality.</p>
HR source ¹⁾	<ul style="list-style-type: none"> – ECG (default) – derives the heart rate from the ECG signal. – Arterial pressure – derives the heart rate from the arterial blood pressure signal. The heart rate parameter field label changes to APR and appears in the color of ART. – SpO₂ – derives the heart rate from the pulse oximetry signal. The heart rate parameter field label changes to PLS and appears in the color of SpO₂. – Auto – derives the heart rate either from the ECG signal or other available sources. If an ECG signal is not available, the M540 switches to Arterial pressure, and then to SpO₂. 	Selects a different source for the heart rate when the ECG channel is unavailable due to artifact resulting from surgical procedures.

Selection	Available settings	Description
Show all leads	None	Shows all ECG waveforms. Press anywhere in the waveform area to access additional ECG waveforms. Press Menu to close all the ECG waveforms.
Size all ECG ¹⁾	0.25, 0.5, 1 (default), 2, 4, 8 mV/cm	Sets the amplitude of ALL displayed ECG leads.
Color ¹⁾	Red, White, Yellow, Green (default), Light blue, Blue, Purple, Orange	Determines the color of the ECG waveforms, and the arrhythmia/ST parameter labels and values.
Settings – ECG 2 page		
Pacer detection (Not available in neonatal mode)	<ul style="list-style-type: none"> – On (default) – Off – the message Pacer off appears in the waveform channel – Fusion – the message Pacer fusion appears in the waveform channel 	Determines whether pacer impulses are detected.
QRS sync marker	<ul style="list-style-type: none"> – On – displays QRS synchronization markers – Off (default) 	Determines whether vertical white markers appear on the waveform to identify QRS complexes. The markers help determine when it is safe to perform synchronized cardioversion.
Cable type ¹⁾ <i>(TruST is only available with a 6-wire lead set)</i>	<ul style="list-style-type: none"> – Auto (default) – 3, 5, 6, and 12 leads (if activated) <p>When using the ECG extension cable, the system always assumes the cable is a 6-wire lead set.</p>	When set to Auto , it detects the number of connected lead wires automatically. If auto mode does not detect the connected lead set, it allows you to select the cable type manually. "12" denotes a combination of a 6-wire lead set and 4-wire lead set for 12-lead monitoring.
ARR lead 1 ¹⁾ ARR lead 2 ¹⁾	ECG1, ECGII (arrhythmia lead 1 default), ECGIII, ECGaVR, ECGaVL, ECGaVF, ECGV (arrhythmia lead 2 default), ECGV+, ECGV1 to ECGV6	Assigns the lead for QRS processing.
ARR processing ¹⁾	ECG1, ECG1&2 (default) The ECG1&2 selection is not available if the neonatal patient category is selected.	ECG1 setting – arrhythmia processing occurs only on the lead selected as arrhythmia lead 1. ECG1&2 setting – arrhythmia processing occurs on the leads selected as arrhythmia lead 1 and arrhythmia lead 2.
¹⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.		
QRS threshold	<ul style="list-style-type: none"> – Normal (default) – Low 	This function is only available for adult and pediatric patients. Normal – detects QRS complexes ≥ 0.5 mV. Low – detects QRS complexes ≥ 0.2 mV.





Selection	Available settings	Description
Rest ECG setup page		
Gender	<ul style="list-style-type: none"> – Unknown (default) – Male – Female 	
Race	<ul style="list-style-type: none"> – Unknown (default) – Caucasian – Asian – African – Other 	

Respiration parameter setup functions

All respiration parameter setup functions take place in the respiration dialog.

Selection	Available settings	Description
Settings 1 page		
Resp. lead ¹⁾	I, II (default)	Selects the lead for respiration monitoring.
Relearn	None	Initiates a relearning of the respiration signal, only in Auto mode.
Marker ¹⁾	On, Off (default)	Superimposes a vertical line on the respiration waveform when a breath is detected.
Monitoring ¹⁾	<ul style="list-style-type: none"> – On (default in neonatal mode) – Off (default in adult/pediatric mode) 	Turns respiration monitoring on or off.
Apnea time ¹⁾	Off , 10, 15 (default), 20, 25, 30 s	Determines how long an apnea has to last before an alarm is triggered.
Apnea archive ¹⁾	<ul style="list-style-type: none"> – Off – Str/Rec – a recording as well as an event storage is triggered automatically in response to an apnea. – Store (default) – a waveform segment is stored in response to an apnea. – Record – a recording is triggered automatically in response to an apnea. 	<p>Determines what happens in response to an apnea.</p> <p>In case of false apnea alarms, it is advised to observe the patient's breathing pattern (belly or chest), and reposition electrodes accordingly, or to adjust the detection threshold manually.</p>
Color ¹⁾	Red, White, Yellow, Green, Light blue (default), Blue, Purple, Orange	Determines the color of the waveforms, and the parameter labels and values.
Change parameter	A list of currently available parameters.	Changes the parameter field to a different parameter.
¹⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.		

Selection	Available settings	Description
Settings 2 page		
Mode¹⁾	Auto (default), Manual	Determines the processing mode for the breath-related impedance change.
Coincidence¹⁾	On, Off (default)	Determines whether you are alerted when the respiratory rate is within 20% of the heart rate, which is an indication that the M540 is counting heart beats as respiration.
Resp. threshold	10% to 100% (in 10% increments) – default: 50%	Adjusts the breath detection threshold.
¹⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.		



SpO2 and Pulse CO-Ox monitoring with Masimo SET MCable

All SpO2 parameter setup functions take place in the SpO2 dialog.

Selection	Available settings	Description
Tone volume ¹⁾	Off , 5, 10 (default), 20, 30, 40, 50, 60, 70, 80, 90, 100%	Adjusts the volume of the pulse tone. If you dock the M540 in an IACS configuration, this setting is replaced by the pulse tone volume setting of the Cockpit. When you undock the M540, this setting is replaced by the Transport pulse tone setting configured under the Alarm setup tab.
Tone source ¹⁾	<ul style="list-style-type: none"> – ECG (default) – the heart blip pulsates with each detected pulse. – SpO2 	Selects the source of the pulse tone which affects the ECG and SpO2 parameter field display. For the SpO2 selection, the higher the pitch of the tone, the higher the SpO2 saturation percentage.
Bar graph ¹⁾	On , Off (default)	Displays a bar graph that is proportional to the pulse rate and strength.
Averaging time ¹⁾	2 to 4, 4 to 6, 8 (default), 10, 12, 14, 16 s	Determines how quickly the reported SpO2 value responds to changes in the patient's oxygen saturation. A longer averaging time provides a more accurate result. However, in clinical situations where rapid physiological changes have to be monitored, use a shorter averaging time
Sensitivity mode ¹⁾	<ul style="list-style-type: none"> – Normal (default) – standard mode – APOD (adaptive probe off detection) – the least sensitive mode for detecting a reading on patients with low perfusion. Provides the best detection for detached sensors. This mode is useful for patients at particular risk for sensors becoming detached such as children or patients who are restless. – Max. – provides maximum sensitivity for poor signals (this mode is recommended for patients with low perfusion or when the low perfusion or low signal quality message is displayed in APOD or Normal sensitivity mode. Max. mode is not recommended for care areas where patients are not monitored visually, such as general wards. It is designed to interpret and display data at the measuring site when the signal may be weak due to decreased perfusion.) 	Determines the level of detection sensitivity.
Fast SAT mode ¹⁾	On , Off (default) When the Averaging time setting is set to 2 to 4 s or 4 to 6 s, the Fast SAT mode selection is grayed out.	Activates rapid tracking of arterial oxygen saturation changes.

Selection	Available settings	Description
Color ¹⁾	Red, White (default), Yellow, Green, Light blue, Blue, Purple, Orange	Determines the color of the waveforms, and the parameter labels and values.
Change parameter	A list of currently available parameters.	Changes the parameter field to a different parameter.
¹⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.		

SpO2 and pulse rate monitoring with Nellcor Oximax MCable

All SpO2 parameter setup functions take place in the SpO2 dialog.

Selection	Available settings	Description
Tone volume ¹⁾	Off , 5, 10 (default), 20, 30, 40, 50, 60, 70, 80, 90, 100%	Sets the volume of the pulse tone. <ul style="list-style-type: none"> – If you dock the M540 in an IACS configuration, this setting is replaced by the pulse tone volume setting of the Cockpit. – When you undock the M540, this setting is replaced by the Transport pulse tone setting configured under the Alarm setup tab
Tone source ¹⁾	– ECG (default) – the heart rate blip pulsates with each detected pulse. – SpO2	Selects the source of the pulse tone which affects the ECG and SpO2 parameter field display. For the SpO2 selection, the higher the pitch of the tone, the higher the SpO2 saturation percentage.
Bar graph ¹⁾	On, Off (default)	Displays a bar graph that is proportional to the pulse rate and strength.
Response mode ¹⁾	– Normal (default) – 90% change within 5 to 7 seconds – Fast – 90% change within 2 to 4 s	Establishes the frequency the oximeter uses to calculate, record, and display SpO2 saturation levels. <ul style="list-style-type: none"> – Normal mode responds to changes in blood oxygen saturation in 5 to 7 seconds – Fast mode responds to changes in blood oxygen saturation levels in 2 to 4 seconds when calculating %SpO2.
SatSeconds ¹⁾	Off (default), 10, 25, 50, 100	This selection does the following: <ul style="list-style-type: none"> – Analyzes desaturation events by multiplying their duration (seconds) by the number of percentage points the patient exceeds the alarm limit. – Eliminates nuisance alarms caused by brief and numerous violations of lower and upper alarm limits. – Overrides the alarm validation setting and the SpO2 high priority desaturation alarm for neonatal patients.
Color ¹⁾	Red, White (default), Yellow, Green, Light blue, Blue, Purple, Orange	Determines the color of the waveforms, and the parameter labels and values.
Change parameter	A list of currently available parameters.	Changes the parameter field to a different parameter.
¹⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.		

Temperature parameter setup functions

All temperature parameter setup functions take place in the temperature dialog.

Selection	Available settings	Description
Ta Tb	<ul style="list-style-type: none"> – TOral – TEso – TNasal – TRect – TBlad – Tcore – TBlid1 – TBlinkt – TSkin – TR – TL 	<p>Ta configures the first temperature value on the M540.</p> <p>Tb configures the second temperature value on the M540.</p>
ΔT		Difference Ta–Tb
T1a T1b	<ul style="list-style-type: none"> – T1Oral – T1Eso – T1Nasal – T1Rect – T1Blad – T1core – T1Blid1 – T1Blinkt – T1Skin – T1R – T1L 	<p>T1a configures the third temperature value.</p> <p>T1b configures the fourth temperature value.</p>
ΔT1		Difference T1a–T1b
Color ¹⁾	Red, White (default), Yellow, Green, Light blue, Blue, Purple, Orange	Determines the color of the parameter labels and values.
Change parameter	A list of currently available parameters.	Changes the parameter field to a different parameter.
¹⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.		

Non-invasive blood pressure parameter setup functions

All non-invasive blood pressure parameter setup functions take place in the non-invasive blood pressure dialog.

Selection	Available settings	Description
Settings		
Interval time ¹⁾ (Cardiac bypass mode automatically deactivates interval measurements)	Off (default), 1 min, 2 min, 2.5 min, 3 min, 5 min, 10 min, 15 min, 20 min, 25 min, 30 min, 45 min, 60 min, 120 min, 240 min	Defines intervals for the non-invasive blood pressure measurements.
Inflation mode	Adult (default), Pediatric , Neonate	Sets a threshold for maximum cuff inflation.
Continuous mode ¹⁾	On , Off (default)	Initiates successive non-invasive blood pressure measurements for 5 min.
Chime ¹⁾	On , Off (default)	Determines if a tone sounds at the end of a completed non-invasive blood pressure measurement.
Venous stasis ¹⁾	On , Off (default)	Stops blood flow to the lower part of the cuffed limb for a fixed time.
Color ¹⁾	Red , White (default), Yellow , Green , Light blue , Blue , Purple , Orange	Determines the color of the parameter labels and values.
Change parameter	A list of currently available parameters.	Changes the parameter field to a different parameter.
¹⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.		

To assign a pressure label manually

- 1 Touch the invasive pressure parameter field.
 - 2 Touch the **Settings** tab.
 - 3 Touch **Edit label**.
 - 4 Touch the appropriate pod label **Pod 1A label**, **Pod 1B label**, and so on.
 - 5 Touch the new pressure label.
 - 6 Touch **X** to close the dialog.
- Or
- 1 Touch the **Menu** function key.
 - 2 Touch **Label IP**.

Zeroing a specific transducer


To zero a specific sensor

- 1 Touch the appropriate invasive pressure parameter field.
- 2 Touch the **Settings** tab.
- 3 Align the transducer to the level of the heart (phlebostatic axis point, fifth intercostal space and midaxillary line).
- 4 Close the transducer stopcock to the patient and open it to air.
- 5 Touch **Zero** on the M540.
If the zeroing of the transducer is successful, the message **Zero accepted** appears. If zeroing fails, the message **did not zero** appears. In that case, repeat steps three to five.
- 6 Touch **X** to close the dialog.

Zeroing all pressure transducers

To zero all pressure transducers from the hemodynamic pods

- 1 Align the transducer to the level of the heart (phlebostatic axis point, fifth intercostal space and midaxillary line).
- 2 Close the stopcocks to the patient, and open them to air.
- 3 Press the $\Rightarrow 0 \Leftarrow$ key (A) on the Hemo4, Hemo2, or the MPod – QuadHemo.



- 4 Verify that the transducers have been zeroed. If zeroing failed, repeat steps two and three.

Invasive pressure parameter setup functions

All invasive pressure setup functions take place in the invasive pressure dialog.

The limits dialog contains the **Auto set** and **Alarm** buttons for configuring the alarm functions.

Selection	Available settings	Description
Settings		
Zero	None	Zeroes only the pressure indicated on the invasive pressure page and displays the time and date of the last zeroing (see page 20).
Edit label ¹⁾	<p>ART, AOR, FEM, AXL, RAD, UAP, BRA, LA, LV, PA, RV, RA, ABD, BDP, CVP, ESO, FEMV, ICP, ICP2, ICP3, ICP4, LA, GPM, RAD, UVP, GP1 to GP8</p> <p>The defaults are as follows:</p> <ul style="list-style-type: none"> – Channel 1: GP1 – Channel 2: GP2 – Channel 3: GP3 – Channel 4: GP4 – Channel 5: GP5 – Channel 6: GP6 – Channel 7: GP7 – Channel 8: GP8 	Allows the user to assign a label to each pressure channel 1 through 8.
Filter ¹⁾	8 and 16 Hz (default)	Selects the filter setting applied to the invasive pressure signal.
Color ¹⁾	<p>Red, White, Yellow, Green, Light blue, Blue, Purple, Orange</p> <p>The various invasive pressure parameters have the following defaults:</p> <ul style="list-style-type: none"> – ART, AOR, FEM, AXL, RAD, UAP, BRA, GP1 to GP4 = Red – PA, LV, BDP = Yellow – CVP, ABD, ESO, FEMV, GPM, UVP = Blue – ICP, ICP2, ICP3, ICP4, LA = Purple – RA, RV, GP5 to GP8 = Orange 	Determines the color of the waveforms, and the parameter labels and values.
Change parameter	A list of currently available parameters.	Changes the parameter field to a different parameter.
<p>¹⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.</p>		

Mainstream CO₂ monitoring

CO₂ parameter setup

All setup functions for CO₂ parameters take place in the CO₂ dialog.

Selection	Available settings	Description
Zero <i>(only available if a CO₂ device is connected)</i>	None	Zeroes the CO ₂ sensor if necessary. The CO ₂ sensor stores a new zero point for CO ₂ measurements.
Atm. pressure ¹⁾	570 to 800 mmHg 760 mmHg (default)	Determines the ambient pressure setting of the sensor and compensates for pressure effects. Failure to compensate for pressure can cause inaccurate measurements.
Gas compens. ¹⁾	Air (default), N₂O/O₂ , O₂ > 50% , HeliOx	Compensates for supplemental oxygen or N ₂ O or HeliOx . Failure to compensate for supplemental oxygen can cause inaccurate measurements.
RRc apnea time ²⁾	Off (default), 10, 15, 20, 25, 30 s	Specifies the time the M540 waits before reporting a cessation of breathing as an apnea event.
Apnea archive ²⁾	Off, Store (default), Str/Rec, Record	Determines what happens in response to an apnea.
Airway adapter	Reusable (default), Disposable	Determines the type of airway adapter used for CO ₂ monitoring. Compensates for the type of airway adapter that is being used. Requires the user to match the adapter with the configuration setting at the M540; if the adapters do not match, the CO ₂ value displayed is compromised.
Color ²⁾	Red, White, Yellow (default), Green, Light blue, Blue, Purple, Orange	Determines the color of the waveforms, and the parameter labels and values.
Change parameter	CO ₂ (default) Examples: HR, SpO ₂ , PLS CO-Ox, CO ₂ , NIBP, RRI, T, T1, GP1, GP2, GP3, GP4, ST	Changes the parameter field to a different parameter.
¹⁾ This setting is a user default that is identical for all patient categories and is also part of the profile. ²⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.		

Microstream CO₂ monitoring

CO₂ parameter setup

Setup functions for CO₂ parameters take place in the CO₂ dialog within the **Microstream** tab.

Selection	Available settings	Description
RRc apnea time ¹⁾	Off (default) 10 s, 15 s, 20 s, 25 s, 30 s	Specifies how long the M540 waits before reporting a cessation in breathing as an apnea event.
Apnea archive ¹⁾	Off, Store (default), STr/Rec	Determines what happens in response to an apnea event.
Next service in:	Informational only (settings are not applicable)	The remaining number of hours until maintenance is required.
Averaging	Instantaneous Last valid breath 10 s, 20 s (default), 30 s	Controls the specific time or the interval used to select the maximum measured etCO₂ and the minimum measured inCO₂ .
Change parameter	CO₂ (default) ECG, ST, NIBP , etc.	Changes the current parameter to another parameter.
Color ¹⁾	Red, White, Yellow (default), Green, Light blue, Blue, Purple, Orange	Determines the color of the CO ₂ waveform, and the parameter labels and values.
Last calibration:	Informational only	Displays the date of the last calibration.
¹⁾ This setting is a patient default which may be unique for each patient category; it is part of the profile.		

This quick reference guide only applies to
Infinity® M540 VG7.n
with the Serial No.:

If no Serial No. has been filled in by Dräger, this quick reference guide is provided for general information only and is not intended for use with any specific machine or unit.

This document is provided for customer information only, and will not be updated or exchanged without customer request.



Directive 93/42/EEC
concerning Medical Devices



The radio equipment in the Infinity M540 patient monitor complies with the Radio Equipment Directive (2014/53/EU). A copy of the Declaration of Conformity is available at the following Internet address: www.draeger.com/doc-radio

 Manufacturer

Dräger Medical Systems, Inc.

3135 Quarry Road
Telford, PA 18969-1042
U.S.A.

(215) 721-5400
(800) 4DRAGER
(800) 437-2437
(215) 723-5935
<http://www.draeger.com>



In Europe, Middle East, Africa, Latin America,
Asia Pacific distributed by

Drägerwerk AG & Co. KGaA

Moislinger Allee 53 - 55
D-23542 Lübeck
Germany
+49 451 8 82-0
+49 451 8 82-20 80
<http://www.draeger.com>

MS34794 – RI 01 en
© Drägerwerk AG & Co. KGaA
Edition: 1 – 2018-09

Dräger reserves the right to make modifications to the equipment without prior notice.

