


AIRVO

Your patient has been commenced on humidified oxygen via the AIRVO system in conjunction with CCO.

The aim of which is to improve respiratory function and patient comfort.

This equipment can deliver both high flows and high percentages of oxygen.

To attain what flow and percentage the patient is receiving please press only the  button.

Please:

- **DO NOT** disconnect oxygen hose from wall.
- **DO NOT** turn off or unplug power
- **DO NOT** adjust flow via oxygen flowmeter. The flowmeter is for oxygen delivery but the actual oxygen concentration of oxygen given is display in the equipment screen. A ward patient should not be on more than 50% oxygen unless discuss with senior doctor.
- **DO NOT** attempt to change settings without CCO support

**FOR ANY ISSUES AT ALL PLEASE
CONTACT CCO VIA 7000**

All disposables excluding irrigation water
can be obtained from GHDU